

Product Datasheet

MED13L antibody (orb167715)

Description	Goat polyclonal antibody to MED13L.
Species/Host	Goat
Reactivity	Bovine, Canine, Human, Porcine
Conjugation	Unconjugated
Tested Applications	ELISA
Target	MED13L (aa361-375)
Preservatives	Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH 7.3 with 0.5% bovine serum albumin. Aliquot and store at -20°C. Minimize freezing and thawing.
Storage	Aliquot and store at -20°C. Minimize freezing and thawing.
Note	For research use only
Application notes	WB: Preliminary experiments gave an approx 80kDa band in lysates of cell line HeLa after 1µg/ml antibody staining. Please note that currently we cannot find an explanation in the literature for the band we observe given the calculated size of 243kDa according
Protein Sequence	SPKRSGKIPPKLHNH
Clonality	Polyclonal
MW	243
Entrez	23389
Dilution Range	ELISA: 1:128000, WB: 1 µg/ml
Expiration Date	12 months from date of receipt.