



Product Datasheet

PRTN3 antibody (orb167717)

Description	Goat polyclonal antibody to PRTN3.
Species/Host	Goat
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, IHC, WB
Target	proteinase 3 / myeloblastin (aa88-98)
Preservatives	Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH 7.3 with 0.5% bovine serum albumin. Aliquot and store at -20°C. Minimize freezing and thawing.
Storage	Aliquot and store at -20°C. Minimize freezing and thawing.
Note	For research use only
Application notes	WB: Approx 28kDa band observed in Human Spleen lysates (calculated MW of 27.8kDa according to NP_002768.3). Recommended concentration: 0.3-1 µg/ml. Primary incubation was 1 hour. IHC: Paraffin embedded Human Spleen. Recommended concentration: 5 µg/ml.
Protein Sequence	HNVRTQEPTQQ
Clonality	Polyclonal
MW	27,8
Entrez	5657
Dilution Range	ELISA: 1:8000, WB: 0.3-1 µg/ml, IHC-P: 5 µg/ml
Expiration Date	12 months from date of receipt.



Western blot analysis of Human Spleen ly...



Immunohistochemical staining of Human Sp...