

Product Datasheet

GlyRS monoclonal antibody (orb1677173)

Description	Mouse IgG1 kappa. Liquid in PBS, pH 7.3, 30% glycerol, and 0.01% sodium azide.
Species/Host	Mouse
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	FC, IF, WB
Immunogen	Recombinant fusion protein of human GlyRS. The exact sequence is proprietary.
Preservatives	Mouse IgG1 kappa. Liquid in PBS, pH 7.3, 30% glycerol, and 0.01% sodium azide.
Concentration	1mg/ml
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze-thaw cycles.
Note	For research use only
Clonality	Monoclonal
MW	~ 85 kDa
Uniprot ID	P41250
Dilution Range	WB (1/1000 - 1/2000), IF/ICC (1/10 - 1/50)
Expiration Date	12 months from date of receipt.