

Product Datasheet

Gamma-catenin monoclonal antibody (orb1677267)

Description	Mouse IgG3. Supplied in crude ascites with 0.01% sodium azide.
Species/Host	Mouse
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	KLH-conjugated synthetic peptide encompassing a sequence within the center region of human Gamma-catenin. The exact sequence is proprietary.
Preservatives	Mouse IgG3. Supplied in crude ascites with 0.01% sodium azide.
Concentration	1mg/ml
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze-thaw cycles.
Note	For research use only
Clonality	Monoclonal
MW	~ 83 kDa
Uniprot ID	P14923
Dilution Range	WB (1/500 - 1/1000)
Expiration Date	12 months from date of receipt.