

Product Datasheet

SIRT6 monoclonal antibody (orb1677496)

Description	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide, pH 7.3.
Species/Host	Mouse
Reactivity	Human, Monkey, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	IF, IP, WB
Immunogen	Purified recombinant human SIRT6 protein expressed in E.coli.
Preservatives	Liquid in PBS containing 50% glycerol, 0.5% rAlbumin and 0.02% sodium azide, pH 7.3.
Concentration	1mg/ml
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze-thaw cycles.
Note	For research use only
Clonality	Monoclonal
MW	~ 42 kDa
Uniprot ID	Q8N6T7
Dilution Range	WB (1/500 - 1/1000), IF/ICC (1/50 - 1/100), IP (1/10 - 1/50)
Expiration Date	12 months from date of receipt.