

## Product Datasheet

### RAD51A monoclonal antibody (orb1677643)

|                            |                                                                                          |
|----------------------------|------------------------------------------------------------------------------------------|
| <b>Description</b>         | Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide, pH 7.3.          |
| <b>Species/Host</b>        | Mouse                                                                                    |
| <b>Reactivity</b>          | Human, Monkey, Mouse, Rat                                                                |
| <b>Conjugation</b>         | Unconjugated                                                                             |
| <b>Tested Applications</b> | WB                                                                                       |
| <b>Immunogen</b>           | Recombinant protein of human Rad51.                                                      |
| <b>Preservatives</b>       | Liquid in PBS containing 50% glycerol, 0.5% rAlbumin and 0.02% sodium azide, pH 7.3.     |
| <b>Concentration</b>       | 1mg/ml                                                                                   |
| <b>Storage</b>             | Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze-thaw cycles. |
| <b>Note</b>                | For research use only                                                                    |
| <b>Clonality</b>           | Monoclonal                                                                               |
| <b>MW</b>                  | ~ 37 kDa                                                                                 |
| <b>Uniprot ID</b>          | <a href="#">Q06609</a>                                                                   |
| <b>Dilution Range</b>      | WB (1/500 - 1/1000)                                                                      |
| <b>Expiration Date</b>     | 12 months from date of receipt.                                                          |