

Product Datasheet

JNK1/2/3 (Phospho-T183/Y185) Rabbit monoclonal antibody (orb1677737)

Description	Liquid in 50mM Tris-Glycine (pH 7.4), 0.15M NaCl, 50% Glycerol, 0.01% Sodium azide and 0.05% BSA.
Species/Host	Rabbit
Reactivity	Human, Rat
Conjugation	Unconjugated
Tested Applications	IP, WB
Immunogen	A synthetic phosphopeptide corresponding to residues surrounding Thr183/Tyr185 of human JNK1
Preservatives	Liquid in 50mM Tris-Glycine (pH 7.4), 0.15M NaCl, 50% Glycerol, 0.01% Sodium azide and 0.05% rAlbumin.
Concentration	1mg/ml
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze-thaw cycles.
Note	For research use only
Clonality	Monoclonal
MW	~ 46,50 kDa
Uniprot ID	P45983
Dilution Range	WB (1/500 - 1/1000), IP (1/10 - 1/50)
Expiration Date	12 months from date of receipt.