

## Product Datasheet

### Filamin A/B (Phospho-S2152/2107) Rabbit monoclonal antibody (orb1677738)

|                            |                                                                                                        |
|----------------------------|--------------------------------------------------------------------------------------------------------|
| <b>Description</b>         | Liquid in 50mM Tris-Glycine (pH 7.4), 0.15M NaCl, 50% Glycerol, 0.01% Sodium azide and 0.05% BSA.      |
| <b>Species/Host</b>        | Rabbit                                                                                                 |
| <b>Reactivity</b>          | Human, Rat                                                                                             |
| <b>Conjugation</b>         | Unconjugated                                                                                           |
| <b>Tested Applications</b> | IF, IHC, WB                                                                                            |
| <b>Immunogen</b>           | A synthetic phosphopeptide corresponding to residues surrounding Ser2152 of human Filamin A            |
| <b>Preservatives</b>       | Liquid in 50mM Tris-Glycine (pH 7.4), 0.15M NaCl, 50% Glycerol, 0.01% Sodium azide and 0.05% rAlbumin. |
| <b>Concentration</b>       | 1mg/ml                                                                                                 |
| <b>Storage</b>             | Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze-thaw cycles.               |
| <b>Note</b>                | For research use only                                                                                  |
| <b>Clonality</b>           | Monoclonal                                                                                             |
| <b>MW</b>                  | ~ 280 kDa                                                                                              |
| <b>Uniprot ID</b>          | <b>O75369</b>                                                                                          |
| <b>Dilution Range</b>      | WB (1/500 - 1/1000), IHC (1/50 - 1/100), IF/ICC (1/50 - 1/100)                                         |
| <b>Expiration Date</b>     | 12 months from date of receipt.                                                                        |