

Product Datasheet

Histone H3 (AcK27) Rabbit monoclonal antibody (orb1677740)

Description

Liquid in 50mM Tris-Glycine (pH

0.01% Sodium azide and 0.05% BSA.

Species/Host

Rabbit

Reactivity

Human

Conjugation

Unconjugated

Tested

ChIP, IF, IHC, IP, WB

Applications
Immunogen

A synthetic acetyl-peptide corresponding to residues surrounding Lys27 of human Histone H3

Preservatives

Liquid in 50mM Tris-Glycine (pH 7.4), 0.15M NaCl, 50% Glycerol, 0.01% Sodium azide and 0.05% rAlbumin.

Concentration

1mg/ml

Storage

Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze-thaw cycles.

Note

For research use only

Clonality

Monoclonal

MW

~ 15 kDa

Uniprot ID
P68431
Dilution Range

WB (1/500 - 1/1000), IHC (1/50 - 1/100), IF/ICC (1/50 - 1/100), IP (1/10 - 1/50), ChIP (1/10 - 1/50)

Expiration Date

12 months from date of receipt.