

Product Datasheet

EIF2S1 (Phospho-S51) Rabbit monoclonal antibody (orb1677853)

Description	Liquid in 50mM Tris-Glycine (pH 7.4), 0.15M NaCl, 50% Glycerol, 0.01% Sodium azide and 0.05% BSA.
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	A synthetic phosphopeptide corresponding to residues surrounding Ser51 of human EIF2S1
Preservatives	Liquid in 50mM Tris-Glycine (pH 7.4), 0.15M NaCl, 50% Glycerol, 0.01% Sodium azide and 0.05% rAlbumin.
Concentration	1mg/ml
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze-thaw cycles.
Note	For research use only
Clonality	Monoclonal
MW	~ 42 kDa
Uniprot ID	P05198
Dilution Range	WB (1/500 - 1/1000)
Expiration Date	12 months from date of receipt.