

Product Datasheet

DNA Polymerase lambda polyclonal antibody (orb1678086)

Description

Liquid in 0.42% Potassium

pH 7.3, 30% glycerol, and 0.01% sodium a...

Species/Host

Rabbit

Reactivity

Human

Conjugation

Unconjugated

Tested

IHC, WB

Applications
Immunogen

KLH-conjugated synthetic peptide encompassing a sequence within the C-term region of human DNA Polymerase lambda. The exact sequence is proprietary.

Preservatives

Liquid in 0.42% Potassium phosphate, 0.87% Sodium chloride, pH 7.3, 30% glycerol, and 0.01% sodium azide.

Concentration

1mg/ml

Storage

Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze-thaw cycles.

Note

For research use only

Clonality

Polyclonal

MW

~ 59 kDa

Uniprot ID

Q9UGP5

Dilution Range

WB (1/500 - 1/1000), IHC (1/50 - 1/200)

Expiration Date

12 months from date of receipt.