

Product Datasheet

Recombinant Human Creatine kinase-MM(CK-MM) antigen (orb1678443)

Description
Recombinant Human cKM-M

expression system with a His tag at the N-term...

Species/Host

Mouse

Tested

CLIA, EIA, LF

Applications
Target

CK-MM

Preservatives

20 mM PB, 150 mM NaCl, 0.1% Proclin 300,pH7.4

Form/Appearance

Supplied as a 0.2 µm filtered solution of 20 mM PB, 150 mM NaCl, 5% glycerol, pH 7.4

Concentration

Lot specific(±20%)

Storage

Short term (up to 7 days) store at 2-8°C. Long term, store at -20°C. Avoid multiple freeze/thaw cycles.

Note

For research use only

Application notes

Use: Calibration/quality control.
Quality control: Mol Mass: 42.8 kDa
AP Mol Mass: 43 kDa Purity: ≥90% as determined by reducing SDS-PAGE

Isotype

IgG1 kappa

Purity

> 95% as determined by SDS-PAGE.

Source

E.coli

Uniprot ID

P05413

Expiration Date

6 months from date of receipt.