

Product Datasheet

Recombinant Human Cardiac Troponin I-C(cTnI+C) antigen (orb1678445)

Description
Recombinant Human Cardiac

produced by E.coli expression system and t...

Species/Host

Mouse

Tested Applications

CLIA, EIA, LF

Target

cTnl+C

Preservatives

20 mM PB, 150 mM NaCl, 0.1% Proclin 300,pH7.4

Form/Appearance

Supplied as a 0.2 µm filtered solution of 20 mM PB, 150 mM NaCl, 5% glycerol, pH 7.4

Concentration

Lot specific(±20%)

Storage

Short term (up to 7 days) store at 2-8°C. Long term, store at -20°C. Avoid multiple freeze/thaw cycles.

Note

For research use only

Application notes

Use: Calibration/quality control.
Quality control: Mol Mass: 62.48 kDa AP Mol Mass: 62 kDa Purity: ≥90% as determined by reducing SDS-PAGE

Isotype

IgG1 kappa

Purity

> 95% as determined by SDS-PAGE.

Source

E.coli

Uniprot ID

P16860

Expiration Date

6 months from date of receipt.