

Product Datasheet

RBP4 Antibody Pair (orb1678461)

Description	RBP4 Antibody Pair
Reactivity	Canine, Equine, Gallus, Human, Rabbit
Tested Applications	ELISA
Preservatives	50% Glycerol, 0.01M PBS, PH 7.4
Form/Appearance	Liquid
Note	For research use only
Application notes	We recommend using the capture antibody at a concentration of 0.6ug/ml and the detection antibody at a concentration of 1ug/ml.Optimal dilutions should be determined experimentally by the researcher.
Isotype	Capture Isotype: IgG1. Detection Isotype: IgG1
Uniprot ID	P02753
Expiration Date	12 months from date of receipt.