

Product Datasheet

PLA2G7 Antibody Pair (orb1678463)

Description	PLA2G7 Antibody Pair
Reactivity	Human
Tested Applications	ELISA
Preservatives	50% Glycerol, 0.01M PBS, PH 7.4
Form/Appearance	Liquid
Note	For research use only
Application notes	We recommend using the capture antibody at a concentration of 1ug/ml and the detection antibody at a concentration of 2ug/ml. Optimal dilutions should be determined experimentally by the researcher.
Isotype	Capture Isotype: IgG1. Detection Isotype: IgG
Uniprot ID	Q13093
Expiration Date	12 months from date of receipt.