

## Product Datasheet

### MT1E Antibody Pair (orb1678503)

|                            |                                                                                                                                                                                                            |
|----------------------------|------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| <b>Description</b>         | MT1E Antibody Pair                                                                                                                                                                                         |
| <b>Reactivity</b>          | Bovine, Canine, Equine, Goat, Human, Porcine, Rabbit                                                                                                                                                       |
| <b>Tested Applications</b> | ELISA                                                                                                                                                                                                      |
| <b>Preservatives</b>       | 50% Glycerol, 0.01M PBS, PH 7.4                                                                                                                                                                            |
| <b>Form/Appearance</b>     | Liquid                                                                                                                                                                                                     |
| <b>Note</b>                | For research use only                                                                                                                                                                                      |
| <b>Application notes</b>   | We recommend using the capture antibody at a concentration of 0.5ug/ml and the detection antibody at a concentration of 0.5ug/ml. Optimal dilutions should be determined experimentally by the researcher. |
| <b>Isotype</b>             | Capture Isotype: IgG. Detection Isotype: IgG                                                                                                                                                               |
| <b>Uniprot ID</b>          | <b>P04732</b>                                                                                                                                                                                              |
| <b>Expiration Date</b>     | 12 months from date of receipt.                                                                                                                                                                            |