

Product Datasheet

Egln3 Antibody Pair (orb1678505)

Description	Egln3 Antibody Pair
Reactivity	Equine, Gallus, Goat, Mouse, Porcine, Rabbit, Rat
Tested Applications	ELISA
Preservatives	50% Glycerol, 0.01M PBS, PH 7.4
Form/Appearance	Liquid
Note	For research use only
Application notes	We recommend using the capture antibody at a concentration of 0.5ug/ml and the detection antibody at a concentration of 1ug/ml. Optimal dilutions should be determined experimentally by the researcher.
Isotype	Capture Isotype: IgG. Detection Isotype: IgG
Uniprot ID	Q91UZ4
Expiration Date	12 months from date of receipt.