

Product Datasheet

ERH Antibody Pair (orb1678527)

Description	ERH Antibody Pair
Reactivity	Gallus, Human
Tested Applications	ELISA
Preservatives	50% Glycerol, 0.01M PBS, PH 7.4
Form/Appearance	Liquid
Note	For research use only
Application notes	We recommend using the capture antibody at a concentration of 0.5ug/ml and the detection antibody at a concentration of 0.2ug/ml. Optimal dilutions should be determined experimentally by the researcher.
Isotype	Capture Isotype: IgG. Detection Isotype: IgG
Uniprot ID	P84090
Expiration Date	12 months from date of receipt.