

Product Datasheet

RPS10 Antibody Pair (orb1678534)

Description	RPS10 Antibody Pair
Reactivity	Canine, Gallus, Guinea pig, Human, Porcine
Tested Applications	ELISA
Preservatives	50% Glycerol, 0.01M PBS, PH 7.4
Form/Appearance	Liquid
Note	For research use only
Application notes	We recommend using the capture antibody at a concentration of 1ug/ml and the detection antibody at a concentration of 0.125ug/ml.Optimal dilutions should be determined experimentally by the researcher.
Isotype	Capture Isotype: IgG. Detection Isotype: IgG
Uniprot ID	P46783
Expiration Date	12 months from date of receipt.