

Product Datasheet

GPR161 Antibody Pair (orb1678555)

Description	GPR161 Antibody Pair
Reactivity	Bovine, Canine, Equine, Goat, Guinea pig, Human, Porcine, Rat, Sheep
Tested Applications	ELISA
Preservatives	50% Glycerol, 0.01M PBS, PH 7.4
Form/Appearance	Liquid
Note	For research use only
Application notes	We recommend using the capture antibody at a concentration of 0.8ug/ml and the detection antibody at a concentration of 0.34ug/ml. Optimal dilutions should be determined experimentally by the researcher.
Isotype	Capture Isotype: IgG. Detection Isotype: IgG
Uniprot ID	Q8N6U8
Expiration Date	12 months from date of receipt.