

Product Datasheet

P4HB Antibody Pair (orb1678568)

Description	P4HB Antibody Pair
Reactivity	Bovine, Canine, Human, Mouse, Porcine, Sheep
Tested Applications	ELISA
Preservatives	50% Glycerol, 0.01M PBS, PH 7.4
Form/Appearance	Liquid
Note	For research use only
Application notes	We recommend using the capture antibody at a concentration of 0.9ug/ml and the detection antibody at a concentration of 0.25ug/ml. Optimal dilutions should be determined experimentally by the researcher.
Isotype	Capture Isotype: IgG. Detection Isotype: IgG
Uniprot ID	P07237
Expiration Date	12 months from date of receipt.