

Product Datasheet

BCAP31 Antibody Pair (orb1678572)

Description	BCAP31 Antibody Pair
Reactivity	Bovine, Canine, Guinea pig, Human, Mouse, Sheep
Tested Applications	ELISA
Preservatives	50% Glycerol, 0.01M PBS, PH 7.4
Form/Appearance	Liquid
Note	For research use only
Application notes	We recommend using the capture antibody at a concentration of 0.2ug/ml and the detection antibody at a concentration of 0.125ug/ml. Optimal dilutions should be determined experimentally by the researcher.
Isotype	Capture Isotype: IgG. Detection Isotype: IgG
Uniprot ID	P51572
Expiration Date	12 months from date of receipt.