

Product Datasheet

PSME2 Antibody Pair (orb1678580)

Description	PSME2 Antibody Pair
Reactivity	Bovine, Canine, Gallus, Goat, Hamster, Human, Mouse, Porcine, Rabbit, Sheep
Tested Applications	ELISA
Preservatives	50% Glycerol, 0.01M PBS, PH 7.4
Form/Appearance	Liquid
Note	For research use only
Application notes	We recommend using the capture antibody at a concentration of 1.6ug/ml and the detection antibody at a concentration of 0.34ug/ml. Optimal dilutions should be determined experimentally by the researcher.
Isotype	Capture Isotype: IgG. Detection Isotype: IgG
Uniprot ID	Q9UL46
Expiration Date	12 months from date of receipt.