

Product Datasheet

EIF4A2 Antibody Pair (orb1678590)

Description	EIF4A2 Antibody Pair
Reactivity	Bovine, Canine, Gallus, Goat, Guinea pig, Human, Mouse, Porcine, Rabbit, Rat, Sheep
Tested Applications	ELISA
Preservatives	50% Glycerol, 0.01M PBS, PH 7.4
Form/Appearance	Liquid
Note	For research use only
Application notes	We recommend using the capture antibody at a concentration of 1.8ug/ml and the detection antibody at a concentration of 0.5ug/ml. Optimal dilutions should be determined experimentally by the researcher.
Isotype	Capture Isotype: IgG. Detection Isotype: IgG
Uniprot ID	Q14240
Expiration Date	12 months from date of receipt.