

Product Datasheet

Recombinant Carcinoembryonic Antigen (CEA) (orb1678791)

Description	Recombinant protein , expressed in Escherichia coli (E. coli) and purified by standard chromatograph...
Tested Applications	SDS-PAGE, WB
Preservatives	20mM Tris, 150mM NaCl, pH 8.0, containing 0.01% SKL, 5% Trehalose.
Storage	Shipped at 2-8°C. Aliqote and Storet at -80°C; Avoid repeated freeze/thaw cycles.
Note	For research use only
Purity	> 90% pure as determined by Silver stained SDS-PAGE analysis.
MW	45KDa
Expiration Date	6 months from date of receipt.