

Product Datasheet

Mouse Anti Estrogen Receptor Alpha Antibody (orb1678887)

Description	The product is affinity purified and specifically recognizes the Human Estrogen Receptor alpha. The ...
Conjugation	Unconjugated
Immunogen	Recombinant Human Estrogen Receptor Alpha
Preservatives	PBS, pH 7.4, containing 0.02% NaN ₃ , 50% glycerol.
Form/Appearance	Liquid
Note	For research use only
Clonality	Monoclonal
Expiration Date	12 months from date of receipt.