

Product Datasheet

Mouse Anti Interleukin 6 (IL6) Antibody (orb1678918)

Description	"The product is affinity purified and specifically recognizes the Human Interleukin 6 (IL6). The a...
Conjugation	Unconjugated
Immunogen	Recombinant IL6 (Val30~Met212) expressed in 293F cell
Preservatives	0.01M PBS, pH 7.4, containing 0.05% Proclin-300, 50% glycerol.
Form/Appearance	Liquid
Note	For research use only
Clonality	Monoclonal
Expiration Date	12 months from date of receipt.