



Product Datasheet

IgE Antibody (APC) (orb1678979)


Description	Immunoglobulin E (IgE) is a 180 kDa soluble protein serving as an antigen-specific unit of mast cell...
Reactivity	Human
Conjugation	APC
Tested Applications	FC
Immunogen	Purified human IgE.
Target	IgE
Preservatives	Stabilizing phosphate buffered saline (PBS), pH 7.4, 15 mM sodium azide
Concentration	0.1 mg/ml
Storage	Store at 2-8°C. Protect from prolonged exposure to light. Do not freeze.
Note	For research use only
Application notes	Flow cytometry: Recommended dilution: 0.5-3 µg/ml
Isotype	Mouse IgG1
Clonality	Monoclonal
Expiration Date	12 months from date of receipt.



Flow
cytometry
multicolor
surface
staini...



Separation of
human IgE
positive
CD45dim...



Flow
cytometry
surface
staining
pattern ...