



Product Datasheet

CD85g Antibody (APC) (orb1678980)

Description	CD85g / ILT7 (immunoglobulin-like transcript 7) is a cell surface protein that is expressed on plasm...
Reactivity	Human
Conjugation	APC
Tested Applications	FC
Target	CD85g
Preservatives	Stabilizing phosphate buffered saline (PBS), pH 7.4, 15 mM sodium azide
Storage	Store at 2-8°C. Protect from prolonged exposure to light. Do not freeze.
Note	For research use only
Application notes	Flow cytometry: The reagent is designed for analysis of human blood cells using 10 µl reagent / 100 µl of whole blood or 10 ⁶ cells in a suspension. The content of a vial (1 ml) is sufficient for 100 tests.
Isotype	Mouse IgG1 kappa
Clonality	Monoclonal
Uniprot ID	P59901
Entrez	23547
Expiration Date	12 months from date of receipt.



Separation of human CD85g positive leuko...



Flow cytometry surface staining pattern ...