

Product Datasheet

alpha 2a Adrenergic Receptor Immunizing peptide (orb16877)

Description	alpha 2a Adrenergic Receptor Immunizing peptide
Preservatives	100µg of Lyophilised peptide, formulation De-ionized (DI) water.
Storage	Shipped at ambient temperature, store at -20°C
Note	For research use only
Application notes	This peptide was used as immunogen to produce orb16877 (alpha 2a Adrenergic Receptor antibody)
Protein Sequence	TERRP NGLGPERS
Expiration Date	12 months from date of receipt.