

## Product Datasheet

### AKR7A3 antibody (orb169579)

|                          |                                              |
|--------------------------|----------------------------------------------|
| <b>Description</b>       | Mouse monoclonal antibody to AKR7A3          |
| <b>Species/Host</b>      | Mouse                                        |
| <b>Reactivity</b>        | Human                                        |
| <b>Conjugation</b>       | Unconjugated                                 |
| <b>Note</b>              | For research use only                        |
| <b>Application notes</b> | Monoclonal Antibodies, Anti-Human<br>Enzymes |
| <b>Clonality</b>         | Monoclonal                                   |
| <b>Expiration Date</b>   | 12 months from date of receipt.              |