

## Product Datasheet

### TAS1R3 Immunizing peptide (orb16993)

|                          |                                                                             |
|--------------------------|-----------------------------------------------------------------------------|
| <b>Description</b>       | TAS1R3 Immunizing peptide                                                   |
| <b>Preservatives</b>     | 100µg of Lyophilised peptide, formulation<br>De-ionized (DI) water.         |
| <b>Storage</b>           | Shipped at ambient temperature, store at -<br>20°C                          |
| <b>Note</b>              | For research use only                                                       |
| <b>Application notes</b> | This peptide was used as immunogen to<br>produce orb16993 (TAS1R3 antibody) |
| <b>Protein Sequence</b>  | GQNDGNTGNQ GKHE                                                             |
| <b>Expiration Date</b>   | 12 months from date of receipt.                                             |