


Product Datasheet

Mouse TNFRSF13B/TACI ELISA Kit (orb170149)

Description	Mouse TNFRSF13B/TACI ELISA Kit
Reactivity	Mouse
Range	46.9 pg/ml - 3,000 pg/ml
Tested Applications	ELISA
Immunogen	Expression system for standard: NS0; Immunogen sequence: A2-C128
Storage	Store at 4°C for 6 months, at -20°C for 12 months. Avoid multiple freeze-thaw cycles (Ships with gel ice, can store for up to 3 days in room temperature. Freeze upon receiving.)
Note	For research use only
Sample Types	cell culture supernatants, cell lysates, serum and plasma (heparin, EDTA).
MW	26947 MW
Uniprot ID	Q9ET35
Sensitivity	10 pg/ml
Dilution Range	Dilution ratio of 1:1, concentration in serum and plasma is less than the lowest standard, 46.9 pg/ml.
Expiration Date	Please enquire.



Mouse
TNFRSF13B/TACI
PicoKine ELISA
Kit ...