


Product Datasheet

Mouse CD163 ELISA Kit (orb170151)

Description	Mouse CD163 ELISA Kit
Reactivity	Mouse
Range	0.78 ng/ml - 50 ng/ml
Tested Applications	ELISA
Immunogen	Expression system for standard: NS0; Immunogen sequence: V39-T1045
Storage	Store at 4°C for 6 months, at -20°C for 12 months. Avoid multiple freeze-thaw cycles (Ships with gel ice, can store for up to 3 days in room temperature. Freeze upon receiving.)
Note	For research use only
Sample Types	cell culture supernatants, serum and plasma (heparin).
MW	120921 MW
Uniprot ID	Q2VLH6
Sensitivity	20 pg/ml
Dilution Range	Dilution ratio of 1:1, concentration in serum and plasma is 6-8 ng/ml.
Expiration Date	Please enquire.



Mouse
CD163
PicoKine
ELISA Kit
standard ...