


Product Datasheet

Mouse TNFSF12/Tweak ELISA Kit (orb170156)

Description	Mouse TNFSF12/Tweak ELISA Kit
Reactivity	Mouse
Range	15.6 pg/ml - 1,000 pg/ml
Tested Applications	ELISA
Immunogen	Expression system for standard: E.coli; Immunogen sequence: S94-H249
Storage	Store at 4°C for 6 months, at -20°C for 12 months. Avoid multiple freeze-thaw cycles (Ships with gel ice, can store for up to 3 days in room temperature. Freeze upon receiving.)
Note	For research use only
Sample Types	cell culture supernatants, cell lysates, serum and plasma (heparin, EDTA).
MW	27441 MW
Uniprot ID	054907
Sensitivity	10 pg/ml
Dilution Range	Dilution ratio of 1:1, concentration in serum and plasma is around 70 pg/ml.
Expiration Date	Please enquire.



Mouse
TNFSF12/TWEAK
PicoKine ELISA
Kit s...