


Product Datasheet

Mouse IL-1 RL1 / ST2 ELISA Kit (orb170159)

Description	Mouse IL-1 RL1 / ST2 ELISA Kit
Reactivity	Mouse
Range	156 pg/ml - 10,000 pg/ml
Tested Applications	ELISA
Immunogen	Expression system for standard: NS0; Immunogen sequence: S27-R332
Storage	Store at 4°C for 6 months, at -20°C for 12 months. Avoid multiple freeze-thaw cycles (Ships with gel ice, can store for up to 3 days in room temperature. Freeze upon receiving.)
Note	For research use only
Sample Types	cell culture supernatants, cell lysates, serum and plasma (heparin, EDTA).
MW	64801 MW
Uniprot ID	P14719
Sensitivity	10 pg/ml
Dilution Range	Dilution ratio of 1:20, concentration in serum and plasma is around 21-77 ng/ml.
Expiration Date	Please enquire.



Mouse
IL1RL1/ST2
PicoKine
ELISA Kit
stan...