


Product Datasheet

Mouse SP-D ELISA Kit (orb170163)

Description	Mouse SP-D ELISA Kit
Reactivity	Mouse
Range	312 pg/ml - 20,000 pg/ml
Tested Applications	ELISA
Immunogen	Expression system for standard: CHO; Immunogen sequence: A20-F374
Storage	Store at 4°C for 6 months, at -20°C for 12 months. Avoid multiple freeze-thaw cycles (Ships with gel ice, can store for up to 3 days in room temperature. Freeze upon receiving.)
Note	For research use only
Sample Types	cell culture supernatants, cell lysates, serum and plasma (heparin, EDTA).
MW	37688 MW
Uniprot ID	P50404
Sensitivity	10 pg/ml
Dilution Range	Dilution ratio of 1:1, concentration in serum and plasma is around 1.5 ng/ml.
Expiration Date	Please enquire.



Mouse SP-D
PicoKine
ELISA Kit
standard c...