


Product Datasheet

Mouse ACE2 ELISA Kit (orb170165)

Description	Mouse ACE2 ELISA Kit
Reactivity	Mouse
Range	62.5 pg/ml - 4,000 pg/ml
Tested Applications	ELISA
Immunogen	Expression system for standard: CHO; Immunogen sequence: Q18-T740
Storage	Store at 4°C for 6 months, at -20°C for 12 months. Avoid multiple freeze-thaw cycles (Ships with gel ice, can store for up to 3 days in room temperature. Freeze upon receiving.)
Note	For research use only
Sample Types	cell culture supernatants, serum and plasma (heparin).
MW	92368 MW
Uniprot ID	Q8R0I0
Sensitivity	10 pg/ml
Dilution Range	Dilution ratio of 1:1, concentration in serum and plasma is 0.6-1 ng/ml.
Expiration Date	Please enquire.



Mouse ACE2
PicoKine
ELISA Kit
standard c...