


Product Datasheet

Mouse Progranulin ELISA Kit (orb170168)

Description	Mouse Progranulin ELISA Kit
Reactivity	Mouse
Range	0.78 ng/ml - 50 ng/ml
Tested Applications	ELISA
Immunogen	Expression system for standard: NS0; Immunogen sequence: T18-L589
Storage	Store at 4°C for 6 months, at -20°C for 12 months. Avoid multiple freeze-thaw cycles (Ships with gel ice, can store for up to 3 days in room temperature. Freeze upon receiving.)
Note	For research use only
Sample Types	cell culture supernatants, cell lysates, serum and plasma (heparin, EDTA).
MW	63458 MW
Uniprot ID	P28798
Sensitivity	10 pg/ml
Dilution Range	Dilution ratio of 1:100, concentration in serum and plasma is 492-2197 ng/ml.
Expiration Date	Please enquire.



Mouse
Progranulin
PicoKine
ELISA Kit
sta...