


Product Datasheet

Mouse Betacellulin/BTC ELISA Kit (orb170169)

Description	Mouse Betacellulin/BTC ELISA Kit
Reactivity	Mouse
Range	7.8 pg/ml - 500 pg/ml
Tested Applications	ELISA
Immunogen	Expression system for standard: E.coli; Immunogen sequence: D32-Y111
Storage	Store at 4°C for 6 months, at -20°C for 12 months. Avoid multiple freeze-thaw cycles (Ships with gel ice, can store for up to 3 days in room temperature. Freeze upon receiving.)
Note	For research use only
Sample Types	cell culture supernatants, cell lysates, serum and plasma (heparin, EDTA).
MW	19664 MW
Uniprot ID	Q05928
Sensitivity	10 pg/ml
Dilution Range	Dilution ratio of 1:1, concentration in serum and plasma is less than the lowest standard, 7.8 pg/ml.
Expiration Date	Please enquire.



Mouse
Betacellulin/BTC
PicoKine ELISA
Ki...