


Product Datasheet

Mouse MFGE8 / Lactadherin / Milk Fat Globule 1 ELISA Kit (orb170170)

Description	Mouse MFGE8 / Lactadherin / Milk Fat Globule 1 ELISA Kit
Reactivity	Mouse
Range	62.5 pg/ml - 4,000 pg/ml
Tested Applications	ELISA
Immunogen	Expression system for standard: NS0; Immunogen sequence: A23-C463
Storage	Store at 4°C for 6 months, at -20°C for 12 months. Avoid multiple freeze-thaw cycles (Ships with gel ice, can store for up to 3 days in room temperature. Freeze upon receiving.)
Note	For research use only
Sample Types	cell culture supernatants, cell lysates, serum and plasma (heparin, EDTA).
MW	51241 MW
Uniprot ID	P21956
Sensitivity	10 pg/ml
Dilution Range	Dilution ratio of 1:10, concentration in serum and plasma is 3.96-10.87 ng/ml.
Expiration Date	Please enquire.



Mouse
MFGE8/Lactadherin
PicoKine ELISA K...