


Product Datasheet

Human MFGE8/Lactadherin ELISA Kit (orb170175)

Description	Human MFGE8/Lactadherin ELISA Kit
Reactivity	Human
Range	62.5 pg/ml - 4,000 pg/ml
Tested Applications	ELISA
Immunogen	Expression system for standard: NS0; Immunogen sequence: L24-C387
Storage	Store at 4°C for 6 months, at -20°C for 12 months. Avoid multiple freeze-thaw cycles (Ships with gel ice, can store for up to 3 days in room temperature. Freeze upon receiving.)
Note	For research use only
Sample Types	cell culture supernatants, serum, plasma (heparin, EDTA) , saliva, urine and human milk.
MW	43123 MW
Uniprot ID	Q08431
Sensitivity	10 pg/ml
Dilution Range	Dilution ratio of 1:1, concentration in serum and plasma is 1712-8887 pg/ml.
Expiration Date	Please enquire.



Human
MFGE8/Lactadherin
PicoKine ELISA K...