


Product Datasheet

Human SPARC ELISA Kit (orb170184)

Description	Human SPARC ELISA Kit
Reactivity	Human
Range	0.78 ng/ml - 50 ng/ml
Tested Applications	ELISA
Immunogen	Expression system for standard: NS0; Immunogen sequence: A18-I303
Storage	Store at 4°C for 6 months, at -20°C for 12 months. Avoid multiple freeze-thaw cycles (Ships with gel ice, can store for up to 3 days in room temperature. Freeze upon receiving.)
Note	For research use only
Sample Types	cell culture supernatants, serum and human milk.
MW	34632 MW
Uniprot ID	P09486
Sensitivity	20 pg/ml
Dilution Range	Dilution ratio of 1:50, concentration in serum is around 115-236 ng/ml.
Expiration Date	Please enquire.



Human
SPARC
PicoKine
ELISA Kit
standard ...