



Product Datasheet

GLI1 Antibody (orb1710827)

Description	GLI1 Antibody
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	IP, WB
Immunogen	between 1056 and 1106 (C-term)
Target	GLI1
Preservatives	Borate Buffered Saline (BBS) pH 8.2 with 0.1% rAlbumin and 0.09% Sodium Azide
Form/Appearance	Whole IgG
Concentration	100 µg/ml
Storage	Shelf life: 1 year from date of receipt. Storage: 2 - 8°C
Note	For research use only
Application notes	Format: Whole IgG
Clonality	Recombinant
Uniprot ID	P08151
Dilution Range	WB: 1:1,000, IP: 20 µl/1 mg lysate
Expiration Date	12 months from date of receipt.



Detection of human GLI1 by western blot....



Detection of human GLI1 by western blot ...