



Product Datasheet

CD86 Antibody (orb1710828)


Description	CD86 Antibody
Species/Host	Rabbit
Reactivity	Mouse
Conjugation	Unconjugated
Tested Applications	FC, ICC, IP, WB
Immunogen	residues 24-245 (ECD)
Target	CD86
Preservatives	Borate Buffered Saline (BBS) pH 8.2 with 0.09% Sodium Azide, rAlbumin-Free
Form/Appearance	Whole IgG
Concentration	1000 µg/ml
Storage	Shelf life: 1 year from date of receipt. Storage: 2 - 8°C
Note	For research use only
Application notes	Format: Whole IgG
Clonality	Recombinant
Uniprot ID	P42082
Dilution Range	WB: 1:1,000, IP: 6 µl/1 mg lysate, ICC: 1:100 to 1:500. Epitope retrieval with citrate buffer pH6.0 is recommended for FFPE cell sections.
Expiration Date	12 months from date of receipt.



Detection of mouse CD86 by western blot....



Detection of mouse CD86 (shaded) in CH27...



Detection of mouse CD86 by western blot ...