

Product Datasheet

CD45, Leucocyte Common Antigen (LCA); Clones PD7/26 & 2B11 (Ready-To-Use) (orb1711979)

Description	CD45, Leucocyte Common Antigen (LCA); Clones PD7/26 & 2B11 (Ready-To-Use)
Species/Host	Mouse
Reactivity	Human
Conjugation	Unconjugated
Immunogen	Isolated neoplastic cells from T-cell lymphoma (2B11); human peripheral blood lymphocytes maintained in T-cell growth factor (PD7/26).
Note	For research use only
Application notes	Ramos, U-698, or GA-10 cells. Tonsil.
Isotype	IgG1, kappa (2B11); IgG1, kappa (PD7/26).
Clonality	Monoclonal
Expiration Date	12 months from date of receipt.