

Product Datasheet

Lambda Light Chain (B-Cell Marker); Clone LcN-2 & ICO-106 (Concentrate) Antibody (orb1712285)

Description	Lambda Light Chain (B-Cell Marker); Clone LcN-2 & ICO-106 (Concentrate) Antibody
Species/Host	Mouse
Reactivity	Human
Conjugation	Unconjugated
Immunogen	Purified human IgG (LcN-2 & ICO-106)
Note	For research use only
Application notes	293T, Raji or hPBL cells. Tonsil or spleen.
Isotype	IgG2a, kappa (LcN-2); IgG1, kappa (ICO-106)
Clonality	Monoclonal
Expiration Date	12 months from date of receipt.