

Product Datasheet

Granulocyte-Colony Stimulating Factor (G-CSF); Clone SPM468 (Concentrate) Antibody (orb1712349)

Description	Granulocyte-Colony Stimulating Factor (G-CSF); Clone SPM468 (Concentrate) Antibody
Species/Host	Mouse
Reactivity	Human, Monkey
Conjugation	Unconjugated
Immunogen	Nuclei from pokeweed mitogen stimulated human peripheral blood lymphocytes.
Note	For research use only
Application notes	HL60 cells, tonsil, or lymph node.
Isotype	IgG1
Clonality	Monoclonal
Expiration Date	12 months from date of receipt.