

Product Datasheet

Cytokeratin 20; Clone EP23 (Ready-To-Use) Antibody (orb1712436)

Description	Cytokeratin 20; Clone EP23 (Ready-To-Use) Antibody
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Immunogen	Rabbits were injected with a synthetic peptide corresponding to residues near the C-terminus in human CK20 protein.
Note	For research use only
Application notes	Colon for normal tissue and colon cancer for abnormal tissue.
Isotype	Rabbit IgG
Clonality	Monoclonal
Expiration Date	12 months from date of receipt.