

Product Datasheet

c-Myc Oncoprotein; Clone 9E10.3 (Concentrate) Antibody (orb1712471)

Description	c-Myc Oncoprotein; Clone 9E10.3 (Concentrate) Antibody
Species/Host	Mouse
Reactivity	Human
Conjugation	Unconjugated
Immunogen	A synthetic peptide, corresponding to aa 408-439 (AEEQKLISEEDLLRKRREQLKHKLEQLRNSCA) from the C-terminus of human c-Myc, coupled to KLH.
Note	For research use only
Application notes	HL-60 cells or breast carcinoma.
Isotype	IgG1, kappa
Clonality	Monoclonal
Expiration Date	12 months from date of receipt.